

Using Stratagene Single Site-Directed Mutagenesis Kit:

Dilutions:

Primers – 1:100dilution of a 500 μ M stock

Primers can be designed using PrimerX

Reaction:

1.5 μ l	5' Primer
1.5 μ l	3' Primer
1.0 μ l	Plasmid [100ng/ μ l]
5.0 μ l	10x PFU Buffer
38 μ l	ddH ₂ O
2.0 μ l	5mM dNTP
1.0 μ l	PFU Polymerase

50 μ l Total

PCR Cycle:

Temperature ($^{\circ}$ C)	Time	Cycles
95	30 seconds	1
95 52 (~10 $^{\circ}$ below T _m of primers) 67	30 seconds 1 minute ~2 minutes/kb ^a	16
67	5-10 minutes	1
4	15-O/N	

Transformation:

Add 1 μ l of Reaction to 50 μ l XL1-Blue cells

Put on ice for 20 minutes

Heat shock at 42 $^{\circ}$ for 45 seconds

Put back on ice for 1 minute

Directly plate the 50 μ l reaction on LB + antibiotic plate

Incubate O/N at 37 $^{\circ}$

^a12 minutes extension for SXL-RRM12 in pGEX-6p